

GO GO™ SUPER PORTABLE



41.6 LBS.
without batteries



Shown with
front storage
bag (standard)

FDA Class II Medical Device¹



UNFOLD



YOUR



WORLD



Quickly folds flat
to a depth of 12"
without having to
remove the seat.



800.800.8586

[pridemobility.com/GoGoSP](https://www.pridemobility.com/GoGoSP)



PRIDE
MOBILITY®

GO GO® SUPER PORTABLE

Shown with standard front storage bag and optional armrests.

SC15M



FDA Class II Medical Device¹

Black Mint True Blue Orange

- 1 FDA Class II Medical Devices are designed to aid individuals with mobility impairments
- 2 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
- 3 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).
- 4 Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- 5 This accessory must be removed prior to using the travel bag.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure.

** Transport of this item on commercial aircraft is subject to federal regulations and airline policy. Please consult your airline prior to flight to ensure your scooter and its 11.7 Ah lithium battery are permitted.



GO GO® TRAVEL BAG BUNDLE*

Conveniently tow your Go Go® Super Portable scooter through airports and parking lots with this durable, protective travel bag bundle.

*Currently only compatible with the Go Go® Super Portable

FEATURES

- Offers a simple folding process without needing to remove the seat.
- Overall depth is reduced to 12", making transport and storage fast and convenient.
- Weighs only **41.6 lbs.** without the battery
- 11.7 Ah **lithium battery is airline-compliant****
- Front and rear lighting for both seeing and being seen
- Integrated USB port for charging mobile devices
- Stands upright when folded and features a handle for easy towing
- Storage compartment under the seat
- Accepts optional armrests
- Front storage bag



ACCESSORIES

- Arm rests
- Cup holder⁵
- Cell phone holder⁵
- Rear view mirrors (single/pair)⁵
- Travel bag bundle

New and exclusive Charge360 features:

- Requires minimal dexterity
- Magnetically connects to the port in any direction
- "Snaps" into place seamlessly and illuminates to indicate a successful connection

SPECIFICATIONS

WEIGHT CAPACITY 300 lbs.

MAXIMUM SPEED² Up to 4.5 mph

DIMENSIONS³ Length: 37.2"
Width: 21.5"

TURNING RADIUS 51.6"

RANGE WITH 11.7 AH BATTERY^{3,4} Up to 10 miles at 220 lbs.
Up to 8.8 miles at 300 lbs.

WARRANTY

FRAME: 3-year limited
ELECTRONICS: 2-year limited
DRIVE MOTORS: 2-year limited
BATTERIES: 2-year limited

STANDS UPRIGHT



pridemobility.com



PRIDE MOBILITY®